

The background of the slide is a photograph of a laboratory. In the foreground, a large yellow robotic arm is visible, with various cables and sensors attached. In the background, another similar yellow robotic arm is positioned over a piece of laboratory equipment. The floor is covered with a grey and white checkered safety mat. The overall scene is brightly lit, typical of a modern research facility.

# High-throughput screening (ToxCast and Tox21 program) and high-content data sources

February 5<sup>th</sup>, 2024

Kristin Eccles, PhD

Research Scientist, Computational Toxicology Research Group

Health Canada

# Outline



1. The Problem
2. What is High- throughput screening and why do we need it
3. The ToxCast/Tox21 Program
  1. Chemical Selection
  2. HTS Analysis/ Assays
  3. Data Processing
    1. Model Fitting, PODs, and Flags
  4. Data Release/ Data Exploration
    1. EPA: CompTox Dashboard
    2. NTP: Integrated Chemical Environment
  5. Data Interpretation
4. Wednesday: HTS Dose-Response Modeling in R

# The Problem

- The “too many chemicals” problem
- A lot of toxicity data (potency, target tissues, etc.) on few chemicals
- Many data poor chemicals
  - Measured in serum, environmental media, in commerce
  - Lack of toxicity information



# The risk assessment process

Step 1: Hazard Identification

Step 2: Hazard Characterization

Step 3: Exposure Assessment



Step 4: Risk Characterization

Step 5: Risk Mitigation



# Traditional methods for hazard characterization

*in vivo*




## The problem:

- Expensive and time consuming
- Animal welfare issues to consider
- Focus on apical endpoints
- Biological relevance to humans



# Toxicology for the 21<sup>st</sup> century

- New Approach Methods:
  - NAMs are broadly defined as any technology, methodology, approach or combination thereof that can be used to replace, reduce or refine animal toxicity testing and allow for more rapid or effective prioritization and assessment of chemicals.
  - E.g., Cell based assays




Nyfeler et al. 2020



# What is high- throughput screening?

- High throughput screening (HTS) is the use of automated equipment to rapidly test thousands to millions of samples for biological activity at varying levels of biological organisation
- *in vitro* approach rapidly provides screening level data on chemical hazard




Vinken et al., 2017

# Tox21 Program

- USA Federal collaboration among EPA, NIH, including National Center for Advancing Translational Sciences (NCATS) and the National Toxicology Program (NTP) at the National Institute of Environmental Health Sciences (NIEHS), and the Food and Drug Administration (FDA).
- Tox21 aims to:
  - Prioritize substances for further in-depth toxicological evaluation
  - Identify mechanisms of action for further investigation (e.g., toxicity-associated and disease-associated pathways)
  - Develop models that better predict how substances will affect biological responses (predictive toxicology)
  - Employ testing methods using human cells (*in vitro* approaches)
  - Reduce time, effort, and costs associated with testing
  - Contribute to the reduction, refinement, and replacement of animals used in toxicity testing


n assays = ~1500  
n chemicals = ~10K




Richard et al. 2020



# Chemical Selection



Richard et al., 2016



Richard et al., 2020

# HTS Assays/ Analysis (non-exclusive list)

Assay	Assay End Point*	Cell Type	Assay Readout
Cell viability Apoptosis	Intracellular ATP, caspase-3/7	Hek293, Jurkat, HepG2, Sh-SY5Y, SK-N-SH, BJ, HUV-EC-C, MRC-5, mesangial, kidney proximal tubules, N2a, H-4-II-E, NIH3T3	Luminescence
Membrane integrity	LDH release, protease release	HEK293, mesangial	Fluorescence, luminescence
Mitochondrial toxicity	Membrane potential	HepG2	Fluorescence
DNA damage	Micronucleus	CHO	Fluorescence
Cytokine	IL-8, TNF- $\alpha$	THP-1	Homogenous time-resolved fluorescence
<a href="#">Nuclear receptor</a>	AR, ER $\alpha$ , FXR, PPAR $\delta$ , PPAR $\gamma$ , RXR, TR $\beta$ , VDR, GR, hPXR, AhR, rPXR, CAR, ERR	Hek293, HeLa, HepG2	$\beta$ -Lactamase reporter, luciferase reporter
Toxicity pathway	AP-1, HIF-1 $\alpha$ , SIE, NF- $\kappa$ B, HSR, ESRE, ARE/Nrf2, CREB, p53, real- time cytotoxicity and viability	ME-180, HeLa, HepG2, Hek293, CHO, HCT-116, HepG2	$\beta$ -Lactamase reporter, luciferase reporter
hERG channel	Thallium influx	U-2OS	Fluorescence
G-protein-coupled receptor signaling	TRHR, TSHR	Hek293	Fluorescence, homogeneous time- resolved fluorescence
Enzyme activity	AChE	SH-SY5Y	Fluorescence
Developmental pathway	Retinol signaling, Hedgehog/Gli, SBE/Smad	C3H10T1/2, 3T3, HEK 293T	Luciferase reporter, fluorescence


- Nomination process
  - Assays are selected based on their biological and toxicological relevance and adaptability to miniaturization and qHTS screening

# Understanding assay naming convention

- TOX21\_CAR\_Agonist
- LTEA\_HepaRG\_CYP1A1\_dn (DMSO)
- LTEA\_HepaRG\_CYP1A1\_up (DMSO)
- CCTE\_Simmons\_AUR\_TPO\_dn (methimazole)
- APR\_HepG2\_MitoticArrest\_24h\_up
- APR\_HepG2\_MitoticArrest\_24h\_dn
- APR\_HepG2\_MitoticArrest\_72h\_dn
- APR\_HepG2\_MitoticArrest\_72h\_up


# Data Processing

- Turning plate read outs into concentration-response curves




High Throughput  
Screening Assays


R: EPA's  
tcpl pipeline








# Mode fitting: Concentration-response modeling



**Fig. 1** The three models utilized by the tcpl package. The constant model and its associated formula for  $\mu_i$  is shown in orange, the Hill model and its associated formula for  $\mu_i$  is shown in blue, and the gain-loss model and its associated formula for  $\mu_i$  is shown in blue

- best model is chosen using AIC


Point of Departure (PODs)



**Fig. 3** The point-of-departure estimates calculated by the tcpl package. The shaded rectangle represents the baseline region,  $0 \pm 3bmad$ . The dark striped line represents the efficacy cutoff. The vertical lines show where the point-of-departure estimates are defined: the activity concentration at baseline (ACB) in red, activity concentration at cutoff (ACC) in yellow, and activity concentration at 50% maximal activity ( $AC_{50}$ ) in blue


$bmad$  = baseline median absolute deviation

# tcplfit2



# Hit call/ Flagging

- Hitc: interpret the continuous hitc into active or inactive
  - Must meet activity thresholds
- Flag – omit
  - QA/QC
  - Removed the least reproducible curve-fits
    - Very low  $AC_{50}$  (below the screened conc range)
    - Low efficacy (within 1.2-fold of the cutoff)
    - Low n of tested concentrations
    - library(tcpl): ?MC6\_Methods(lvl = 6)



```
ASSAY: AEID1 (TOX21_ERa_BLA_Agonist_ratio_gain)

NAME: Phenylalanine
CHID: 23463 CASRN: 150-30-1
SPID(S): Tox21_110011
M4ID: 3

HILL MODEL (in red):
tp ga gw
val: 0.375 0.915 0.81
sd: 0.119 0.411 0.294

GAIN-LOSS MODEL (in blue):
tp ga gw la lw
val: 0.375 0.915 0.81 3.1 5.19
sd: 0.119 0.412 0.295 670 2340

CNST HILL GNLS
AIC: -18.94 -49.21 -45.21
PROB: 0 0.88 0.12
RMSE: 0.13 0.04 0.04

MAX_MEAN: 0.313 MAX_MED: 0.313 BMAD: 0.0546
COFF: 0.328 HIT-CALL: 0 FITC: 21 ACTP: 1
```

The hill\_tp and gnls\_tp parameters are equal and greater than coff; however, the maximum median value (max\_med) is not greater than the cutoff making the series inactive.

[https://cran.r-project.org/web/packages/tcpl/vignettes/Data\\_processing-Archive\\_tcpl\\_v2.html](https://cran.r-project.org/web/packages/tcpl/vignettes/Data_processing-Archive_tcpl_v2.html)



# Data Release/ Data Exploration

- **Tox21 Toolbox**


- The Tox21 toolbox contains useful tools for accessing and visualizing the Tox21 quantitative high throughput screening (qHTS) 10K library data, as well as integrating with other publicly available data.
- Integrated Chemical Environment (ICE): <https://ice.ntp.niehs.nih.gov>
  - Has extra curation step

- **EPA CompTox Chemicals Dashboard**

- The EPA Computational Toxicology (CompTox) Chemicals Dashboard is a one-stop-shop for chemistry, toxicity and exposure information for over 760,000 chemicals. Data and models within the Dashboard also help with efforts to identify chemicals of most need of further testing and reducing the use of animals in chemical testing. The CompTox Chemicals Dashboard integrates data from various sources including:
- ToxCast and Tox21 - High-throughput chemical screening data
- General Use – Assay Information: <https://comptox.epa.gov/dashboard/>
- ToxCast: <https://www.epa.gov/comptox-tools/exploring-toxcast-data>

# Data Interpretation/ Utilization

- Translating *in vitro* data into real world meaning
  - Using in pathway analysis (e.g., Estrogen and Androgen Receptors)
  - Input into *in silico* models
    - QSAR/ physiochemical properties for read-across
  - *in vitro* to *in vivo* extrapolation (IVIVE) (Dr. Marc Beal's Lecture)



# Future Directions for Tox21 (Phase III: 2014–present)

- **Toxicodynamic Variability in Developmental Neurotoxicity**
  - *Goal: Incorporate genetic variation into cell-based test systems to better understand potential population differences in response to chemicals that may cause toxic neurological effects.*
- **Retrofitting Existing Tox21 High-Throughput Screening Assays with Metabolic Capability**
  - *Goal: To add xenobiotic metabolism capability to existing Tox21 assays so they provide more accurate and informative data regarding in vivo activity.*
- **Cell Line Selection for High Throughput Transcriptomics (HTT)**
  - *Goal: Develop a strategy for selecting maximally diverse cell types/lines to maximally cover biological targets and pathways for high-throughput chemical screening using gene expression (i.e., transcriptomics).*
- **In Vitro Chemical Disposition**
  - *Goal: Understand the impact of chemical disposition within in vitro test systems across a broad range of chemical categories and develop a computational model to predict differences between the “nominal” concentration of a chemical compared with “true” concentration in the media and cells.*

# Wednesday Hands-On Activity

Concentration-response modeling in R with ToxCast/ Tox21 Data

1. Install R from <https://cran.rstudio.com>
2. Install RStudio from:  
<https://www.rstudio.com/products/rstudio/download/#download>
3. Download (or clone) the repository with the materials for  
Wednesday: <https://github.com/kristineccles/BIM4103-HTS-Dose-Response>
  1. R scripts
  2. Data



# References

- Chang, X., Tan, Y. M., Allen, D. G., Bell, S., Brown, P. C., Browning, L., ... & Mumtaz, M. (2022). IVIVE: Facilitating the use of in vitro toxicity data in risk assessment and decision making. *Toxics*, 10(5), 232.
- Filer, D. L., Kothiya, P., Setzer, R. W., Judson, R. S., & Martin, M. T. (2017). tcpl: the ToxCast pipeline for high-throughput screening data. *Bioinformatics*, 33(4), 618-620.
- Richard, A. M., Judson, R. S., Houck, K. A., Grulke, C. M., Volarath, P., Thillainadarajah, I., ... & Thomas, R. S. (2016). ToxCast chemical landscape: paving the road to 21st century toxicology. *Chemical research in toxicology*, 29(8), 1225-1251.
- Richard, A. M., Huang, R., Waidyanatha, S., Shinn, P., Collins, B. J., Thillainadarajah, I., ... & Tice, R. R. (2020). The Tox21 10K compound library: collaborative chemistry advancing toxicology. *Chemical Research in Toxicology*, 34(2), 189-216.
- Vinken, M., Knapen, D., Vergauwen, L., Hengstler, J. G., Angrish, M., & Whelan, M. (2017). Adverse outcome pathways: a concise introduction for toxicologists. *Archives of toxicology*, 91, 3697-3707.